

Flying with your Tandem Insulin Pump



TANDEM[®]
DIABETES CARE

If you plan to travel by air with your insulin pump, some advance planning may make your security screening go smoother. Your device is safe for use during air travel and complies with FAA wireless transmission standards. It is also designed to withstand common electromagnetic interference and can be safely carried through metal detectors.

However, your Tandem insulin pump should **NOT** be put through machines that use X-rays, including airline luggage X-ray machines and full-body scanners. We recommend asking the security agent for an alternative screening method and request to go through a standard metal detector wearing your pump.


It may be helpful to provide a printout of the information below to the security agent with your pump during screening. If you have any questions about travelling with your pump, please call NZMS Diabetes Customer Care at 0508 634 103. We are here for you 24 hours a day, 7 days a week.



Metal Detectors, hand held body scanners



x-rays, including full body scanners and luggage screening machines



Your Choice. Your Life. Your Freedom.[®]

24/7 Technical Support 0508 634 103  
0800 500 226 | nzms@nzms.co.nz | nzmsdiabetes.co.nz

