

t:slim X2™ Insulin Pump

T:LOCK™ INFUSION SETS AND CARTRIDGE

AutoSoft 90

- 90° insertion angle
- Integrated inserter
- Teflon cannula



AutoSoft 30

- 30° insertion angle
- Integrated inserter
- Teflon cannula



TruSteel Set

- Manual straight in insertion
- Stainless steel/metal cannula



VariSoft™ Set

- Manual angled set
- Teflon cannula



t:lock Infusion Set and Cartridge Pharmac Order Codes (10 per pack)

AutoSoft 90 Infusion Set t:lock 6mm/60cm	2556804	TruSteel Infusion Set t:lock 6mm/60cm	2556863
AutoSoft 90 Infusion Set t:lock 6mm/110cm	2556812	TruSteel Infusion Set t:lock 6mm/80cm	2616491
AutoSoft 90 Infusion Set t:lock 9mm/60cm	2556820	TruSteel Infusion Set t:lock 8mm/60cm	2556898
AutoSoft 90 Infusion Set t:lock 9mm/110cm	2556839	TruSteel Infusion Set t:lock 8mm/80cm	2612550
AutoSoft 30 Infusion Set t:lock 13mm/60cm	2556847	VariSoft Infusion Set t:lock 13mm/60cm	2679833
AutoSoft 30 Infusion Set t:lock 13mm/110cm	2556855		

t:lock Cartridge

t:lock Cartridge	2556790
------------------	---------

0800 500 226
diabetes@nzms.co.nz

nzmsdiabetes.co.nz



a Dexcom company

Always read the label and follow the directions for use. Read the warnings available on nzmsdiabetes.co.nz/resources before purchasing. Consult your healthcare professional to see which product is right for you. The t:lock cartridge is an insulin reservoir for Tandem insulin pump. AutoSoft 90, AutoSoft 30, TruSteel infusion sets deliver insulin from Tandem insulin pump, subcutaneously self-inserted. The VariSoft infusion set is for the subcutaneous infusion of insulin. © 2024 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, t:lock, and t:slim X2 are registered trademarks or trademarks of Tandem Diabetes Care, Inc. Distributed in New Zealand by NZMS (2A Fisher Crescent Mt. Wellington, Auckland 1060 New Zealand). NZMS is a Dexcom company. NZDI_TAN_103 October 2024