

Infusion Set Choices

The types of infusion set used on the t:slim X2™ are the same as used on the Animas Vibe however they have different names and some line lengths vary.

We have provided a table below to help you easily identify the changes.

Animas Name	Animas Options	Tandem Name	Tandem Options
Contact Detach	6mm cannula 60cm line 6mm cannula <u>110cm</u> line 8mm cannula 60cm line 8mm cannula <u>110cm</u> line	TruSteel™	6mm cannula 60cm line 6mm cannula <u>81cm</u> line 8mm cannula 60cm line 8mm cannula <u>81cm</u> line
Inset 30	13mm cannula 60cm line 13mm cannula 110cm line	AutoSoft™ 30	13mm cannula 60cm line 13mm cannula 110cm line
Inset II	6mm cannula 60cm line 6mm cannula 110cm line 9mm cannula 60cm line 9mm cannula 110cm line	AutoSoft™ 90	6mm cannula 60cm line 6mm cannula 110cm line 9mm cannula 60cm line 9mm cannula 110cm line